

Examination of signal parameters

Sojung Son

sonsojun@fel.cvut.cz

Examination of signal parameters

From loaded signal/data (voltage, sound, blood pressure, wind speed, ...) examine its basic parameters as amplitude, DC component, peak-peak voltage, rms voltage, frequency content, overshoot, ... Examined parameters depends on type of the signal. Source of the signal is up to you.

→ ECG signal (Electrocardiography)

How does it work?

run the program

search ECG data (xlsx file)

push button to start examination

go through some correction of data

show signal graph and several parameters

parameters : R-R interval, P-R interval, QRS width

end program or reset to examine another data

```

%% Build GUI
hFig = uifigure('Name' , 'Examination of Signal Parameters', ...
    'Position' , [(screenSize(3:4)-figSize)/2 figSize]);

but_searchfile = uibutton(hFig, 'Text', 'Search File', ...
    'Position', [20 530 200 50]);

lab_file = uilabel(hFig, 'Position', [240 530 660 50], ...
    'Text', '', ...
    'BackgroundColor', ones(1,3), ...
    'FontSize', 14);

but_examine = uibutton(hFig, 'Text', 'go',...
    'Position', [920 530 60 50],...
    'Enable', 'off');

axx_signal = uiaxes(hFig, 'Position', [30 90 450 420]);
disableDefaultInteractivity(axx_signal);
axx_signal.Toolbar.Visible = 'off';
lab_property = uilabel(hFig, 'Position', [530 90 450 420], ...
    'Text', '', ...
    'BackgroundColor', ones(1,3), ...
    'FontSize', 16);

sli_xl = uislider(hFig, 'Position', [ 150 70 300 3], ...
    'Orientation', 'horizontal', ...
    'Enable', 'off', ...
    'MajorTickLabels', {}, ...
    'MajorTicks', [],...
    'MinorTicksMode', 'manual', ...
    'Limits', [0 log(10)]);

sli_xp = uislider(hFig, 'Position', [150 50 300 3], ...
    'Orientation', 'horizontal', ...
    'Enable', 'off', ...
    'MajorTickLabels', {}, ...
    'MajorTicks', [], ...
    'MinorTicksMode', 'manual', ...
    'Limits', [-1 1] );

but_restart = uibutton(hFig, 'Text', 'Restart', ...
    'Position', [560 20 200 50], ...
    'Enable', 'off');

but_close = uibutton(hFig, 'Text', 'Close', ...
    'Position', [780 20 200 50]);


lab_sli_xl_1 = uilabel(hFig, 'Text', '-', 'Position', [133 65 20 20]);
lab_sli_xl_2 = uilabel(hFig, 'Text', '+', 'Position', [455 65 20 20]);
lab_sli_xp_1 = uilabel(hFig, 'Text', '<', 'Position', [133 43 20 20]);
lab_sli_xp_2 = uilabel(hFig, 'Text', '>', 'Position', [455 43 20 20]);
lab_x = uilabel(hFig, 'Text', 'X axis', 'Position', [80 60 100 20]);

```


Examination of Signal Parameters

Search File

go



X axis



Restart

Close

%% Appdata

```
setappdata(axx_signal, 'default', 0);
```

```
setappdata(axx_signal, 'new', 0);
```

%% Callbacks

```
but_searchfile.ButtonPushedFcn = @(src,event)choosefile ;
```

```
but_examine.ButtonPushedFcn = @(src,event)startexamine ;
```

```
sli_xl.ValueChangingFcn = @(src,event) xaxis_l(event) ;
```

```
sli_xp.ValueChangingFcn = @(src,event) xaxis_p(event) ;
```

```
but_restart.ButtonPushedFcn = @(src,event)resetall;
```

```
but_close.ButtonPushedFcn = @(src,event)closewindow(hFig) ;
```

%% constants

```
filename = '';
```

```
signal_data = [];
```

```
figSize = [1000, 600];
```


```
screenSize = get(groot, 'ScreenSize');
```

Examination of Signal Parameters

Search File

go

```
function choosefile
    filename = uigetfile('*.xlsx', 'Select Excel File');
    if filename == 0
    else
        set(lab_file, 'Text', filename);
        set(but_examine, 'Enable', 'on');
        signal_data = xlsread(filename);
    end
end
```



X axis

-

<

+

>

Restart

Close

Examination of Signal Parameters

Search File


ECG.xlsx

go

```
function startexamine
    s = size(signal_data);
    if s(1) ~= 1 && s(2) ~= 1
        uialert(hFig, 'Choose valid file.', 'Invalid File');
    else
        %axes
        def_xl = length(signal_data); %default lenght of data
        plot(axx_signal, 1:def_xl, signal_data);
        set(sli_xl, 'Enable', 'on');
        setappdata(axx_signal, 'default', def_xl);
        setappdata(axx_signal, 'XLim', [0 def_xl]);
        setappdata(axx_signal, 'XLimMode', 'manual');
        %label
        [RR_interval, PR_interval, QRS_width] = examine(signal_data);
        text_RR = sprintf('R-R Interval : %d\n normal rate : 60-100 bpm\n\n', RR_interval);
        text_PR = sprintf('P-R Interval : %d\n normally 120-200 ms\n\n', PR_interval);
        text_QRS = sprintf('QRS Width : %d\n normally about 0.12 s\n', QRS_width);
        set(lab_property, 'Text', [text_RR text_PR text_QRS]); %print properties
        %etc
        set([but_searchfile, but_examine], 'Enable', 'off');
        set(but_restart, 'Enable', 'on');
    end
end
```

Close


```
function [RR_interval, PR_interval, QRS_width] = examine(signal_data)
M = max(signal_data);
m = min(signal_data);
thr = (9*M + m)/10;
locs_p = [];
locs_q = [];
locs_r = [];
locs_s = [];
```




```

% R
value_h= 0; % if value > thr
j = 1;
new_i = [];
for i = 2: length(signal_data)-1
    if signal_data(i) <= thr
        if value_h == 1
            new_i(j, 2) = i;
            j = j+1;
            value_h = 0;
        end
    else
        if value_h == 0
            new_i(j, 1) = i;
            value_h = 1;
        end
    end
end
for i = 1:j-1
    rm = max(signal_data(new_i(i, 1):new_i(i, 2)));
    r = find(signal_data(new_i(i, 1):new_i(i, 2)) == rm, 1);
    locs_r = [locs_r, r+new_i(i,1)-1];
end
% R-R interval
RR_interval = fix(mean(diff(locs_r)));

```




```

L = length(locs_r);
j = 1;
for j = 1: L-1
    %S
    sm = min(signal_data(locs_r(j) : locs_r(j)+fix(RR_interval/2)));
    s = find(signal_data(locs_r(j) : locs_r(j)+fix(RR_interval/2)) == sm, 1);
    locs_s = [locs_s, s+locs_r(j)-1];

    %Q
    qm = min(signal_data(locs_r(j) - fix(RR_interval/3) : locs_r(j)));
    q = find(signal_data(locs_r(j) - fix(RR_interval/3) : locs_r(j)) == qm, 1);
    locs_q = [locs_q, q+locs_r(j) - fix(RR_interval/3)-1];
    j = j+1;
end

    %P
    for j = 1 : length(locs_q) % should i change to 2?
        pm = max(signal_data(locs_q(j)-fix(RR_interval/3):locs_q(j)));
        p = find(signal_data(locs_q(j)-fix(RR_interval/3):locs_q(j)) == pm, 1);
        locs_p = [locs_p, p+locs_q(j)-fix(RR_interval/3)-1];
    end
end

```



```

% P-R interval
PR = NaN(1, L-1);
if locs_r(1) > locs_p(1)
    for i = 1:L-1
        PR(i) = locs_r(i) - locs_p(i);
    end
else
    for i = 1: L -2
        PR(i) = locs_r(i+1) - locs_p(i);
    end
end
PR_interval = fix(mean(PR));
% QRS Width
QRS = NaN(1, L-1);
for i = 1:L-1
    QRS(i) = locs_s(i) -locs_q(i);
end
QRS_width = fix(mean(QRS));

```

```

if isempty(locs_r)
    uialert(hFig, 'Choose ECG signal.', 'Invalid Signal Data');
    RR_interval = 'invalid';
    PR_interval = 'invalid';
    QRS_width = 'invalid';
end

```

Examination of Signal Parameters

Search File

ECG.xlsx

go

```
function startexamine
    s = size(signal_data);
    if s(1) ~= 1 && s(2) ~= 1
        uialert(hFig, 'Choose valid file.', 'Invalid File');
    else
        %axes
        def_xl = length(signal_data); %default lenght of data
        plot(axx_signal, 1:def_xl, signal_data);
        set(sli_xl, 'Enable', 'on');
        setappdata(axx_signal, 'default', def_xl);
        setappdata(axx_signal, 'XLim', [0 def_xl]);
        setappdata(axx_signal, 'XLimMode', 'manual');
        %label
        [RR_interval, PR_interval, QRS_width] = examine(signal_data);
        text_RR = sprintf('R-R Interval : %d\n normal rate : 60-100 bpm\n\n', RR_interval);
        text_PR = sprintf('P-R Interval : %d\n normally 120-200 ms\n\n', PR_interval);
        text_QRS = sprintf('QRS Width : %d\n normally about 0.12 s\n', QRS_width);
        set(lab_property, 'Text', [text_RR text_PR text_QRS]); %print properties
        %etc
        set([but_searchfile, but_examine], 'Enable', 'off');
        set(but_restart, 'Enable', 'on');
    end
end
```

Close

